

UREMIC TOXINS

Edited by TOSHIMITSU NIWA

WILEY

UREMIC TOXINS

Edited By

TOSHIMITSU NIWA, MD, PhD Nagoya University Graduate School of Medicine Nagoya, Japan



own-regulated in

hanism of organ calcium. *Miner*

bone metabolism
S201, 2003.
Work Group:
prevention, and
MBD). *Kidney*

3C, Delmar EG, athyroid calcium

i MG, Hercz G, M, Evenepoel P, e TB, Goodman g hemodialysis.

<u>16</u>

β₂-MICROGLOBULIN

Suguru Yamamoto, Junichiro James Kazama, Hiroki Maruyama, Ichiei Narita, and Fumitake Gejyo

16.1 CHEMICAL STRUCTURE AND MOLECULAR WEIGHT

 $β_2$ -microglobulin is a polypeptide of 99 residues that has a molecular weight of 11.8 kDa. It forms the beta chain of the human leukocyte antigen (HLA) class I molecule and has a well-known β-sandwich structure that involves a seven-strand β-pleated structure stabilized with a single disulphide bond (Cys25–Cys80) (Fig. 16.1a). $β_2$ -microglobulin changes conformation under various *in vivo* or *in vitro* conditions. Far-UV spectra show that the fractions of β-sheet and β-turn decrease and the fractions of α-helix and random structure increase with several kinds of treatment, such as acidic pH,¹ 2,2,2-trifluoroethanol,²sodium dodecyl sulfate (SDS),³ lysophospholipids,⁴ nonesterified fatty acids,⁵ heating,⁶ and agitation.⁷ A recent finding shows that the $β_2$ -microglobulin mutant, K58P-W60G, where a Pro residue has been introduced in the type I β-turn, improves chemical and temperature stability and makes folding faster relative to native $β_2$ -microglobulin.⁸ These conformational changes and stability are pivotal for $β_2$ -microglobulin-related amyloid fibril formation/extension (Fig. 16.1b), which induces dialysis-related amyloidosis (DRA).

16.2 METABOLISM AND BIOLOGY

 β_2 -microglobulin is a component of MHC class I molecules, which are present on all nucleated cells. Multiple myeloma, lymphatic leukemia, and malignant lymphoma,

Uremic Toxins, Edited by Toshimitsu Niwa.

^{© 2012} John Wiley & Sons, Inc. Published 2012 by John Wiley & Sons, Inc.



FIGURE 16.1 (a) Structure of β_2 -microglobulin (PDB 1JNJ³⁴). (b) Electron microscopy images of β_2 -miroglobulin-related amyloid fibrils from carpal tunnel tissue in hemodialysis patients. Bar shows 1 μ m.

as well as diseases with activation of the cellular immune system induce increased production of β_2 -microglobulin in serum.

Most β_2 -microglobulin is normally eliminated by the kidney via glomerular filtration and subsequent tubular catabolism. Megalin is a multiligand endocytic receptor involved in the endocytosis into lysosomes and metabolism of a wide variety of glomerular-filtered proteins, including β_2 -microglobulin in the proximal tubule.⁹ These proteins are metabolized in proximal tubular epithelial cells and the metabolites are recovered from general circulation via peritubular capillaries. Thus, proximal tubular injury increases β_2 -microglobulin excretion in urine, while severe kidney damage induces the retention of β_2 -microglobulin in serum due to impaired excretion from the kidney.

16.3 QUANTIFICATION METHODS

Serum level of β_2 -microglobulin is usually measured using immunoassays, such as an indirect solid phase enzyme-linked immunosorbent assay (ELISA). β_2 -microglobulin is the precursor protein for β_2 -microglobulin-related amyloid fibrils and needs conformational change for *in vitro* fibril formation, as described above. Recent clinical studies have attempted to identify the conformational intermediate form of β_2 -microglobulin. Capillary electrophoresis reveals that patients with chronic kidney disease (CKD) undergoing hemodialysis (HD), but not healthy control subjects, have the conformational intermediate form of β_2 -microglobulin in serum.¹⁰ The level of predialysis serum β_2 -microglobulin intermediate was 2.7 ± 1.4 mg/L and native β_2 -microglobulin was 29.4 ± 6.8 mg/L in 31 HD patients. HD using a polymethylme membrane devintermediate f microglobulin treatment. It rr ellular space w tissues (cartila matography-m of lysine-58-c which was fou subjects.¹¹ Ho containing β_2 formational int of onset/progr relationship.

16.4 PLASN SUBJECTS

Advanced CKI impaired metal levels of β_2 -mi to those in norr that the impair HD patients; I microglobulin the precursor continuous hig of carpal tunne

16.5 TOXIC

16.5.1 Adeq

 β_2 -microglobu serum β_2 -micr dialysis for n microglobulin membrane, an

16.5.2 Mort:

In the random β_2 -microglobi



on microscopy 1 hemodialysis

ice increased

a glomerular nd endocytic a wide variety timal tubule.⁹ I the metabo-Thus, proxiwhile severe e to impaired

says, such as λ). β_2 -microid fibrils and bove. Recent diate form of nronic kidney trol subjects, serum.¹⁰ The 1.4 mg/L and HD using a polymethylmethacrylate and online hemodiafiltration (HDF) with a polysulfone (PS) membrane decreased the level of the native form, while any change in the intermediate form was variable.¹⁰ These results indicate that intermediate β_2 -microglobulin is increased in HD patients and is difficult to remove with dialysis treatment. It may suggest that the intermediate form is immobilized in the extracellular space where β_2 -microglobulin-related amyloid has a marked affinity for joint tissues (cartilage, capsule, and synovium). In addition, immunoaffinity-liquid chromatography-mass spectrometry analysis and immunoassay revealed the generation of lysine-58-cleaved and truncated β_2 -microglobulin (Δ K58- β_2 -microglobulin), which was found in serum from 20–40% HD patients but not in serum from control subjects.¹¹ However, this truncated form has not been demonstrated in the tissue containing β_2 -microglobulin-related amyloid.¹² It is not certain whether the conformational intermediate or the truncated form of β_2 -microglobulin has a critical role of onset/progress of DRA, and future studies will be needed to understand the relationship.

16.4 PLASMA/SERUM LEVELS IN UREMIC PATIENTS AND HEALTHY SUBJECTS

Advanced CKD induces to elevate the serum level of β_2 -microglobulin due to the impaired metabolism and excretion in the kidney. The average serum concentration levels of β_2 -microglobulin in patients undergoing HD is significantly higher compared to those in normal subjects (25–45 mg/L vs. 1–2 mg/L).^{10,13–16} It is clearly understood that the impairment of metabolism in the kidney is the main cause of fluid retention in HD patients; however, it is not clearly understood whether the production of β_2 -microglobulin is increased with CKD and/or dialysis treatment. β_2 -microglobulin is the precursor protein of β_2 -microglobulin-related amyloid fibrils in DRA, and continuous higher serum levels of β_2 -microglobulin probably indicate the presence of carpal tunnel syndrome, one of the major symptoms induced by DRA.¹³

16.5 TOXICITY AND CLINICAL RELEVANCE

16.5.1 Adequacy of Dialysis Treatment

 β_2 -microglobulin, a middle-size protein, is a marker for adequacy of dialysis. A high serum β_2 -microglobulin level and a low β_2 -microglobulin *Kt/V* indicate that adequate dialysis for middle-size proteins has not been achieved. Generally, serum β_2 -microglobulin is increased by reduced residual renal function, use of low-flux dialyzer membrane, and long-term dialysis treatment even with high-flux membrane.¹⁵

16.5.2 Mortality

In the randomized Hemodialysis (HEMO) Study, the relationship between serum β_2 -microglobulin levels or dialyzer β_2 -microglobulin clearance and mortality over

252 β₂-MICROGLOBULIN

a period of 2.84 years was analyzed.¹⁵ In time-dependent Cox regression models, predialysis serum B2-microglobulin levels were associated with all-cause mortality after adjustment for residual kidney urea clearance and the number of prestudy years on dialysis.¹⁵ In addition, Okuno et al. reported that higher serum β_2 microglobulin levels showed higher all-cause mortality and noncardiovascular mortality during 3.33 years.¹⁶ In this study, there was no correlation between serum B2-microglobulin and cardiovascular mortality. However, another crosssectional study revealed that serum β_2 -microglobulin is related with heart valve calcification, which is associated with carotid intima-media thickness in HD patients.¹⁴ Furthermore, serum β_2 -microglobulin levels were positively correlated with several cardiovascular risk factors, such as highly-sensitive C-reactive protein, troponin-T, myeloperoxidase, N-terminal pro-B-type natriuretic peptide and inversely correlated with prealbumin and ankle-brachial index.¹⁷ These reports indicate that increased clearance for middle-size molecules, such as β_2 microglobulin, may induce less all-cause motality or less cardiovascular event/ motality.

16.5.3 Dialysis-Related Amyloidosis

DRA is a form of general amyloidosis characterized by the deposition of β_2 microglobulin-related amyloid fibrils.¹⁸ Long-term dialysis patients often have DRA-related clinical manifestations, such as carpal tunnel syndrome and destructive arthropathy associated with cystic bone lesions (Fig. 16.2). The tissue retention of β_2 -microglobulin with dialysis treatment due to the persistence of high serum concentrations for long durations appears to be a prerequisite for the onset of DRA. In addition to the retention of β_2 -microglobulin, partial unfolding of β_2 -microglobulin is believed to be needed for its assembly into B2-microglobulin-related amyloid fibrils, in combination with interaction of disease-related molecules, such as apolipoprotein E,19 glycosaminoglycans,2,20 proteoglycans,²⁰ and phospholipids.⁴ Thus, accumulated β_2 -microglobulin in dialysis patients forms amyloid fibrils with conformational change and stabilization of amyloid fibrils with several disease-related molecules. In addition to deposition of B₂-microglobulin-related amyloid fibrils, inflammation after infiltration of inflammatory cells induces various DRA-related manifestations. In the case of DRA-osteopathy, chronic inflammation is found in the lesion. Inflammatory cells often infiltrate into the synovial and/or disc tissue that contain β_2 microglobulin-related amyloid deposits,²¹ and these inflammatory cells release inflammatory cytokines.²² Thus, amyloid deposition induces local osteolysis through synovial inflammation and subsequent osteoclastogenesis and/or osteoclast activation through three possible pathways: (i) indirect action of inflammatory cytokines through the expression in osteoblasts of receptor activator of nuclear factor-kB ligand/osteoprotegerin ligand (RANKL/OPGL), (ii) direct action of inflammatory cytokines, and (iii) RANKL/OPGL expression in inflammatory cells (Fig. 16.3).²³



FIGURE 16.2 undergone hen showed destruc (b). Congo-red

16.6 THER

16.6.1 Dialy

Therapeutic c dialysis moda relationship is level.^{24,25} Thi microglobulii while the prasignificantly c and dialysis c renal transpla improving Dl

16.6.2 HD

The use of hig microglobulii THERAPEUTIC METHODS TO REMOVE THE TOXINS 253



osition of B2its often have vndrome and g. 16.2). The he persistence rerequisite for bulin, partial issembly into interaction of noglycans,^{2,20} roglobulin in e and stabili-In addition to on after infiltations. In the on. Inflammait contain β_2 -' cells release cal osteolysis and/or osteoof inflammar activator of .), (ii) direct ion in inflam-



FIGURE 16.2 Destructive arthropathy in hip joint due to DRA. A 68 year-old man who had undergone hemodialysis for 34 years had rt. hip joint pain and gait disturbance. Bone X-ray showed destructive arthropathy of rt. femur neck (a), and the man underwent joint replacement (b). Congo-red (c) stained the tissue with orange-green birefringence under polarized light (d).

16.6 THERAPEUTIC METHODS TO REMOVE THE TOXINS

16.6.1 Dialysis Modality

Therapeutic options to improve serum β_2 -microglobulin levels are to use several dialysis modalities, such as HD, HDF, and peritoneal dialysis. A significant inverse relationship is observed between residual renal function and serum β_2 -microglobulin level.^{24,25} This suggests that peritoneal dialysis may keep lower serum levels of β_2 -microglobulin because of better maintenance of intrinsic renal function than HD, while the prevalence of histological DRA in peritoneal dialysis patients is not significantly different from that observed in a group of HD patients matched for age and dialysis duration.²⁶ A radical approach to reducing serum β_2 -microglobulin is renal transplantation, which decreases it from 39 ± 6 to 3.8 ± 1 mg/L as well as improving DRA symptoms and inhibiting the progress of DRA.²⁷

16.6.2 HD with High-Flux Membrane

The use of high-flux dialyzer membrane leads to a reduction in the serum level of β_2 -microglobulin as compared to using low-flux dialyzer membrane. In the HEMO

524 b²-MICKOGLOBULIN



FIGURE 16.3 Possible pathways of bone resorption in DRA-osteopathy. (A) Inflammatory cytokines released from synovial inflitating cells around amyloid deposits act on osteoblasts to express receptor activator of RANKL/OPGL on the cell surface. RANKL/OPGL promote osteoclastic bone resorption through direct contact to osteoclastic lineages. (B) Inflammatory cytokines released from synovial inflitating cells around amyloid deposition act directly on osteoclastic lineages to promote bone resorption without the involvement of RANKL/OPGL on the cell surface, and osteoclastic bone resorption is promoted through direct contact to osteoclastic lineages. (C) Synovial inflitating cells around amyloid depositis express RANKL/OPGL on the cell surface, and osteoclastic bone resorption is promoted through direct contact to osteoclastic lineages. Local administration of osteoprotegerin is expected to prevent osteoclastic bone resorption for through pathways A and C.

Study, the predialysis serum β_2 -microglobulin level was lower in the high-flux membrane group than in the low-flux membrane group (35.3 ± 11.2 mg/L vs. 41.7 ± 11.8 mg/L).¹⁵ In addition, the rate of increase during 36 months was low in the high-flux membrane group.¹⁵ In another study, switching of dialyzer from conventional to high-flux membrane reduced the predialysis serum β_2 -microglobulin level by 22%, from 39.1 ± 11.2 to 30.5 ± 7.2 mg/L.²⁸ The reason why high-flux membrane produces a lower level of serum β_2 -microglobulin is not only that it membrane produces a lower level of serum β_2 -microglobulin is not only that it promotes better clearance, but that it also increases the binding of β_2 -microglobulin to not cells, such as granulocytes, lymphocytes, and monocytes. These cells have to blood cells, such as granulocytes, lymphocytes, and monocytes. These cells have more β_2 -microglobulin on the cell surface in the patients undergoing HD using high-flux more β_2 -microglobulin on the cell surface in the patients undergoing HD using high-flux more β_2 -microglobulin on the cell surface in the patients undergoing HD using high-flux more β_2 -microglobulin on the cell surface in the patients undergoing HD using high-

flux membran significant difi groups.²⁹ Thu: better clearan fewer onsets membrane coi

M CH E.8.81

A β₂-microgle eliminate seru: β₂-microglobulin dependent int molecule. At series before i series before i polysulfone di going HD as (srudy also sho'



flux membrane compared to patients using low-flux membrane, while there is no significant difference in β_2 -microglobulin mRNA expression in the cells among the groups.²⁹ Thus, high-flux membrane decreases serum β_2 -microglobulin levels due to better clearance and less release from the surface of cells. These results indicate fewer onsets of DRA ²⁸ as well as better survival³⁰ with HD using high-flux membrane compared to HD using low-flux membrane.

16.6.3 HD with β_2 -Microglobulin Adsorption Column

A β_2 -microglobulin adsorption column has been developed as a way to directly eliminate serum β_2 -microglobulin, although limitations are encountered in removing β_2 -microglobulin by the improvement of high-flux membrane alone. This adsorption column system is designed for direct hemoperfusion (Fig. 16.4). Adsorption of β_2 -microglobulin by this column is a result both of hydrophobic and molecular size-dependent interactions between the ligand in the column and β_2 -microglobulin molecule. At each standard HD treatment, the adsorption column is connected in series before arterial blood enters the dialyzer. According to a prospective multicenter study, a β_2 -microglobulin adsorption column that was placed in series with a polysulfone dialyzer increased serum β_2 -microglobulin reduction in patients undergoing HD as compared to the control group (74.1 ± 6.1% vs. 60.1 ± 6.3%).³¹ This study also showed improvement of DRA-related symptoms, such as joint pain, and it



FIGURE 16.4 Schematic (a) and photographic (b) representation of blood flow through the dialyzer and β_2 -microglobulin adsorption column.

) Inflammatory on osteoblasts
)PGL promote
) Inflammatory act directly on ANKL/OPGL.
 GL on the cell to osteoclastic eoclastic bone

the high-flux 1.2 mg/L vs. iths was low lialyzer from icroglobulin hy high-flux only that it icroglobulin se cells have D using high-

256 β₂-MICROGLOBULIN

may suggest that the column absorbs not only β_2 -microglobulin, but also other molecules related to inflammation.

16.6.4 Hemodiafiltration

HDF has better clearance of middle-size molecules than HD and is known to reduce the risk of progression of DRA.³² A recent multicenter prospective randomized study examined the treatments, A: HD with low flux for 6 months after high-efficiency online HDF for 6 months or B: high-efficiency online HDF for 6 months.³³ Online HDF showed greater efficiency than HD with low-flux membrane in reducing the basal level of β_2 -microglobulin (22.2 ± 7.8 mg/L vs. 33.5 ± 11.8 mg/L) and of small sized molecules.³³

REFERENCES

- Yamaguchi I, Hasegawa K, Naiki H, Mitsu T, Matuo Y and Gejyo F: Extension of A beta2M amyloid fibrils with recombinant human β₂-microglobulin. *Amyloid* 8:30–40, 2001.
- 2. Yamamoto S, Yamaguchi I, Hasegawa K, Tsutsumi S, Goto Y, Gejyo F, Naiki H: Glycosaminoglycans enhance the trifluoroethanol-induced extension of β_2 -microglobulin-related amyloid fibrils at a neutral pH. *J Am Soc Nephrol* 15:126–133, 2004.
- Yamamoto S, Hasegawa K, Yamaguchi I, Tsutsumi S, Kardos J, Goto Y, Gejyo F, Naiki H: Low concentrations of sodium dodecyl sulfate induce the extension of β₂-microglobulinrelated amyloid fibrils at a neutral pH. *Biochemistry* 43:11075–11082, 2004.
- Ookoshi T, Hasegawa K, Ohhashi Y, Kimura H, Takahashi N, Yoshida H, Miyazaki R, Goto Y, Naiki H: Lysophospholipids induce the nucleation and extension of β₂-microglobulin-related amyloid fibrils at a neutral pH. *Nephrol Dial Transplant* 23:3247–3255, 2008.
- Hasegawa K, Tsutsumi-Yasuhara S, Ookoshi T, Ohhashi Y, Kimura H, Takahashi N, Yoshida H, Miyazaki R, Goto Y, Naiki H: Growth of β₂-microglobulin-related amyloid fibrils by non-esterified fatty acids at a neutral pH. *Biochem J* 416:307–315, 2008.
- Sasahara K, Naiki H and Goto Y: Exothermic effects observed upon heating of β₂microglobulin monomers in the presence of amyloid seeds. *Biochemistry* 45:8760–8769, 2006.
- 7. Sasahara K, Yagi H, Sakai M, Naiki H and Goto Y: Amyloid nucleation triggered by agitation of β_2 -microglobulin under acidic and neutral pH conditions. *Biochemistry* 47:2650–2660, 2008.
- Colombo M, Ricagno S, Barbiroli A, Santambrogio C, Giorgetti S, Raimondi S, Bonomi F, Grandori R, Bellotti V, Bolognesi M: The effects of an ideal β-turn on β₂ microglobulin fold stability. *J Biochem* 150:39–47, 2011.
- 9. Orlando RA, Rader K, Authier F, Yamazaki H, Posner BI, Bergeron JJ, Farquhar MG: Megalin is an endocytic receptor for insulin. J Am Soc Nephrol 9:1759–1766, 1998.
- 10. Uji Y, Motomiya Y and Ando Y: A circulating β_2 -microglobulin intermediate in hemodialysis patients. *Nephron Clin Pract* 111:c173–c181, 2009.

- Corlin DB, S of cleaved β Clin Chem 5
- Giorgetti S, Marini S, Øs 58-cleaved £ fibrils of two 2007.
- Gejyo F, Hor form of amy 314:585–58(
- Ikee R, Hor Association drome, β₂-m Ther Apher
- Cheung AK, Beck GJ, Cl Serum β₂ mi study. J Am
- Okuno S, Is Nishizawa ' maintenance
- 17. Kuragano T, T: The impa in hemodial
- 18. Gejyo F, Yaı M, Hirasaw associated v Biophys Res
- 19. Yamaguchi the depolyn *chemistry* 4
- 20. Yamaguchi N, Yamamo depolymeri: 2003.
- Gejyo F, Mε in β₂-micro
- 22. Ohashi K, F most suscer Int 41:1646
- 23. Kazama JJ, dialysis-rela
- 24. McCarthy J dialysis pat 1994.

but also other

town to reduce domized study high-efficiency onths after HD than HD with 2.2 ± 7.8 mg/L

Extension of A nyloid 8:30-40,

yo F, Naiki H: β_2 -microglobu-33, 2004.

ejyo F, Naiki H: -microglobulin-2004.

H, Miyazaki R, on of β_2 -micro-23:3247–3255,

l, Takahashi N, related amyloid -315, 2008.

heating of β_2 -45:8760–8769,

on triggered by Biochemistry

ondi S, Bonomi 2 microglobulin

Farquhar MG: -1766, 1998. ntermediate in

- Corlin DB, Sen JW, Ladefoged S, Lund GB, Nissen MH and Heegaard NH: Quantification of cleaved β₂-microglobulin in serum from patients undergoing chronic hemodialysis. *Clin Chem* 51:1177–1184, 2005.
- Giorgetti S, Stoppini M, Tennent GA, Relini A, Marchese L, Raimondi S, Monti M, Marini S, Østergaard O, Heegaard NH, Pucci P, Esposito G, Merlini G, Bellotti V: Lysine 58-cleaved β₂-microglobulin is not detectable by 2D electrophoresis in *ex vivo* amyloid fibrils of two patients affected by dialysis-related amyloidosis. *Protein Sci* 16:343–349, 2007.
- Gejyo F, Homma N, Suzuki Y and Arakawa M: Serum levels of β₂-microglobulin as a new form of amyloid protein in patients undergoing long-term hemodialysis. N Engl J Med 314:585–586, 1986.
- Ikee R, Honda K, Oka M, Maesato K, Mano T, Moriya H, Ohtake T, Kobayashi S: Association of heart valve calcification with malnutrition-inflammation complex syndrome, β₂-microglobulin, and carotid intima media thickness in patients on hemodialysis. *Ther Apher Dial* 12:464–468, 2008.
- Cheung AK, Rocco MV, Yan G, Leypoldt JK, Levin NW, Greene T, Agodoa L, Bailey J, Beck GJ, Clark W, Levey AS, Ornt DB, Schulman G, Schwab S, Teehan B, Eknoyan G: Serum β₂ microglobulin levels predict mortality in dialysis patients: Results of the HEMO study. J Am Soc Nephrol 17:546–555, 2006.
- Okuno S, Ishimura E, Kohno K, Fujino-Katoh Y, Maeno Y, Yamakawa T, Inaba M, Nishizawa Y: Serum β₂-microglobulin level is a significant predictor of mortality in maintenance haemodialysis patients. *Nephrol Dial Transplant* 24:571–577, 2009.
- Kuragano T, Kida A, Furuta M, Nanami M, Otaki Y, Hasuike Y, Nonoguchi H, Nakanishi T: The impact of β₂-microglobulin clearance on the risk factors of cardiovascular disease in hemodialysis patients. *ASAIO J* 56:326–332, 2010.
- Gejyo F, Yamada T, Odani S, Nakagawa Y, Arakawa M, Kunitomo T, Kataoka H, Suzuki M, Hirasawa Y, SHirahama T, Cohen AS, Schmid K: A new form of amyloid protein associated with chronic hemodialysis was identified as β₂-microglobulin. *Biochem Biophys Res Commun* 129:701–706, 1985.
- 19. Yamaguchi I, Hasegawa K, Takahashi N, Gejyo F and Naiki H: Apolipoprotein E inhibits the depolymerization of β_2 -microglobulin-related amyloid fibrils at a neutral pH. *Biochemistry* 40:8499–8507, 2001.
- Yamaguchi I, Suda H, Tsuzuike N, Seto K, Seki M, Yamaguchi Y, Hasegawa K, Takahashi N, Yamamoto S, Gejyo F, Naiki H: Glycosaminoglycan and proteoglycan inhibit the depolymerization of β₂-microglobulin amyloid fibrils *in vitro*. *Kidney Int* 64:1080–1088, 2003.
- Gejyo F, Maruyama H, Teramura T, Kazama J, Ei I and Arakawa M: Role of macrophages in β₂-microglobulin-related dialysis amyloidosis. *Contrib Nephrol* 112:97–104, 1995.
- Ohashi K, Hara M, Kawai R, Ogura Y, Honda K, Nihei H, Mimura N: Cervical discs are most susceptible to β₂-microglobulin amyloid deposition in the vertebral column. *Kidney Int* 41:1646–1652, 1992.
- Kazama JJ, Maruyama H and Gejyo F: Osteoclastogenesis and osteoclast activation in dialysis-related amyloid osteopathy. Am J Kidney Dis 38:S156–S160, 2001.
- McCarthy JT, Williams AW and Johnson WJ: Serum β₂-microglobulin concentration in dialysis patients: Importance of intrinsic renal function. J Lab Clin Med 123:495–505, 1994.

258 β₂-MICROGLOBULIN

- Yamamoto S, Kasai A and Shimada H: High peritoneal clearance of small molecules did not provide low serum β₂-microglobulin concentrations in peritoneal dialysis patients. *Perit Dial Int* 23 (Suppl 2): S34–S36, 2003.
- Jadoul M, Garbar C, Vanholder R, Sennesael J, Michel C, Robert A, Noël H, van Ypersele de Strihou C: Prevalence of histological β₂-microglobulin amyloidosis in CAPD patients compared with hemodialysis patients. *Kidney Int* 54:956–959, 1998.
- 27. Mourad G and Argiles A: Renal transplantation relieves the symptoms but does not reverse β_2 -microglobulin amyloidosis. *J Am Soc Nephrol* 7:798–804, 1996.
- Koda Y, Nishi S, Miyazaki S, Haginoshita S, Sakurabayashi T, Suzuki M, Sakai S, Yuasa Y, Hirasawa Y, Nishi T: Switch from conventional to high-flux membrane reduces the risk of carpal tunnel syndrome and mortality of hemodialysis patients. *Kidney Int* 52:1096–1101, 1997.
- 29. Traut M, Haufe CC, Eismann U, Deppisch RM, Stein G and Wolf G: Increased binding of β_2 -microglobulin to blood cells in dialysis patients treated with high-flux dialyzers compared with low-flux membranes contributed to reduced β_2 -microglobulin concentrations. *Results of a cross-over study. Blood Purif* 25:432–440, 2007.
- Locatelli F, Martin-Malo A, Hannedouche T, Loureiro A, Papadimitriou M, Wizemann V, Jacobson SH, Czekalski S, Ronco C, Vanholder R;Membrane Permeability Outcome (MPO) Study Group: Effect of membrane permeability on survival of hemodialysis patients. J Am Soc Nephrol 20:645–654, 2009.
- 31. Gejyo F, Kawaguchi Y, Hara S, Nakazawa R, Azuma N, Ogawa H, Koda Y, Suzuki M, Kaneda H, Kishimoto H, Oda M, Ei K, Miyazaki R, Maruyama H, Arakawa M, Hara M: Arresting dialysis-related amyloidosis: A prospective multicenter controlled trial of direct hemoperfusion with a β₂-microglobulin adsorption column. *Artif Organs* 28:371–380, 2004.
- 32. Nakai S, Iseki K, Tabei K, Kubo K, Masakane I, Fushimi K, Kikuchi K, Shinzato T, Sanaka T, Akiba T: Outcomes of hemodiafiltration based on Japanese dialysis patient registry. *Am J Kidney Dis* 38:S212–S216, 2001.
- 33. Pedrini LA, De Cristofaro V, Comelli M, Casino FG, Prencipe M, Baroni A, Campolo G, Manzoni C, Colì L, Ruggiero P, Acquistapace I, Auriemma L: Long-term effects of high-efficiency on-line haemodiafiltration on uraemic toxicity. A multicentre prospective randomized study. *Nephrol Dial Transplant* 26:2617–2624.
- 34. Verdone G, Corazza A, Viglino P, Pettirossi F, Giorgetti S, Mangione P, Andreola A, Stoppini M, Bellotti V, Esposito G: The solution structure of human β_2 -microglobulin reveals the prodromes of its amyloid transition. *Protein Sci* 11:487–499, 2002.

<u>17</u>

CYTO

BJÖRN ANDI

17.1 CHEN

There is no c cytokines are system, whic response. The (adipocytokin functions. T-I terized by co cytokines and the cytokines middle molec some extent a

An extens possible with as uremic tox Toxin Work C them have by biological me with the acute stabilize IL-6 five different 1 28.0 kDa. The and phosphor

Uremic Toxins, E © 2012 John Wil